

O2 Automatic Change Over Valve

SPECIFICATION SHEET

Features

O2 Automatic Change Over unit provides seamless transition from one cylinder to another without the need to change the cylinders manually, therefore allowing the user to concentrate on the task in hand without the worry of checking individual cylinder contents.

- 0-10 Bar Contents Gauge
- A choice of hose tail, NIST or DISS connections
- Each unit can be custom made to suit individual position and outlet requirements
- Input Pressure 4 Bar
- Output Pressure 4 Bar
- Cylinder change over when pressure drops below 3 Bar

