


Modular Change Over Valve

SPECIFICATION SHEET

Features

The O2 Manual Change Over unit provides seamless transition from one cylinder to another without the need to change the cylinders manually. When the bottle pressure is low on one cylinder simply switch the lever over to use the other cylinder.

- 0-10 Bar Contents Gauge
- A choice of hose tail, NIST or DISS connections
- Each unit can be custom made to suit individual position and outlet requirements
- Input Pressure 4 Bar
- Output Pressure 4 Bar



Description	Part No
O2 – With Hose Tails	1130400V
O2/N2O – With Hose Tails	1130300/2V
O2 – With DISS Connectors	1130401D
O2 – With NIST Connectors	1130401N